

WHILLOCK TARGETED ORAL CARE SYSTEM

The World's Most Effective Way to Keep Your Mouth Healthy



E. B. Whillock, DDS, FICD | (605) 280-3961

WHILLOCK TARGETED ORAL CARE SYSTEM

The World's Most Effective Way to Keep Your Mouth Healthy



People need an easier, more effective method of preventing tooth decay and gum diseases. The Provent system is a simple, proven way to kill and remove disease-causing bacteria that reside between teeth and under the edges of the gums. Conventional toothbrushes and dental floss cannot reach these areas. These revolutionary benefits are achieved using two specially designed interchangeable brush heads: an antimicrobial brush head and a polishing brush head.



Antimicrobial Brush Head

- Electric oral hygiene brush has a rotary-motor drive system that causes adjacent bristle tufts to twist in opposite directions as the bristle tuft action reciprocates back and forth. The unique design of the bristle tufts allows them to deliver antimicrobial agents into areas between teeth, and into the areas between gums and teeth. Additionally, the drive system features a 2½-minute timer, which provides the user with guidance to be more thorough while using the antimicrobial brush head.
- Brush head consists of eight bristle tufts (four on each side) in two different heights. Filaments in the long bristle tufts are slightly stiffer than conventional brushes so they can better enter crevices between teeth. Filaments in short bristle tufts are softer to gently clean crevices along the sides of the teeth.
- Antimicrobial liguigel incorporates a germ-killing agent that when used with the brush kills and removes bacteria living in these critical areas of the mouth. Gel-like liquid is easy to retain in the mouth during brushing. Toothpastes should not be used because they block the rotating tuft mechanisms over a short period of use. The liguigel ensures longevity of the brush head.
- Blue color coordinates with blue colored antimicrobial liguigel.

WHILLOCK TARGETED ORAL CARE SYSTEM

The World's Most Effective Way to Keep Your Mouth Healthy

Polishing Brush Head

- Interchangeable with the same handle; operates using the same rotary-motor drive system.
- Brush head consists of two bristle arrays rotating in opposite directions to polish more efficiently.
- This head is white in color. Coordinates with white polishing/whitening cream. Distinguishes it from the blue Antimicrobial Brush Head.
- Special polishing/whitening cream is used in conjunction with the Polishing Brush Head to help remove stains from the teeth.
- The brush head and the special polishing/whitening cream are especially gentle and nonabrasive to the teeth and gums.



WHILLOCK TARGETED ORAL CARE SYSTEM

The World's Most Effective Way to Keep Your Mouth Healthy





PRODUCT NAME:

WHILLOCK TARGETED ORAL CARE SYSTEM

How to Use the Antimicrobial Brush Head

Step 1: Put about a ½ oz. of the liquigel (about one tablespoon full) into your mouth and keep it in your mouth while you use the brush.

Step 2: Turn the brush on. Place the brush in your mouth. Keep your lips sealed around the shank of the brush to prevent the liquigel from leaking out of your mouth.

Step 3: While your lips are sealed around the shank of the brush, push the tips of the bristles into the spaces between the teeth and along the edge of the gums. Tilt your head to the side you are brushing to pool the area with the germ-killing liquigel. Remember, you are not trying to brush the teeth surfaces with this procedure. You are killing the disease-causing bacteria living in the spaces between the teeth and gums. Do not move the brush back and forth like an ordinary toothbrush. Push the tips of the rotating bristles into the areas between the teeth. Do this about six times or so in one area then move about the width of one tooth to the next area. Work your way around your mouth in small steps while pushing the bristle tips into the areas between the teeth and along the edge of the gums.

Be sure you push the bristles into the spaces along the cheek side and tongue side surfaces of the teeth.

Step 4: After you are sure you have not skipped over any area around the teeth, brush back and forth, and side to side, on the tongue surface. Spit out the liquigel, and rinse the bristles and the head of the brush under running water. There is no need to rinse your mouth; just spit out the remaining liquigel. The entire once-a-day procedure should take less than 2 ½ minutes.

Step 5: Place the brush head back on the handle and store the brush in the charging base.

Keep the charging base plugged in so the brush is fully charged each time you use it.

How to Use the Polishing Brush Head

The polishing brush head & whitening cream are especially gentle and non-abrasive to your teeth & gums. Use the polishing brush head with the polishing/whitening cream to control tooth stains. You can do the germ removal procedure before going to bed at night, or whenever you have time during the day. The cosmetic brushing with the polishing brush head and polishing/whitening cream can be done in the morning or at some other time.



- 
- Kills and removes bacteria from areas under the edges of the gums and between the teeth where other oral-hygiene devices cannot reach.
 - A more effective and efficient way to maintain optimal mouth health at home compared to conventional methods of brushing and flossing. Uses a unique electric brush-head design that works in conjunction with an antimicrobial liquigel.
 - Keeps the mouth healthy and teeth polished: Helps prevent tooth decay and gum disease, and helps decrease the risk of heart disease, stroke and other medical conditions associated with periodontal disease.

- Reduces the time necessary to keep mouth disease-free. (Flossing takes extra time, does not kill bacteria, nor does it reach the bacteria living in the majority of the crevices.)
- Its housing is made of a durable plastic; the contoured handle provides a secure grip during use.

